

## PROSTHETIC REHABILITATION UTILIZING THE HYBRID RESTORATION PLANNING PROTOCOLS



### ***WHAT RECORDS DO WE NEED TO START A CASE?***

#### **WHAT RECORDS DO WE NEED TO START CASE?**

- Photos from side and front at mouth level
  - Lips apart at rest
  - Smile
  - Maximum animation
- Maxillary and mandibular alginate impression / model
- Bite registration
- Face - bow (if desired) SAM III preferred by Lab



### ***WHAT HAPPENS ON THE DAY OF SURGERY?***

#### **WHAT HAPPENS ON THE DAY OF SURGERY?**

- Surgeon extracts teeth, removes bone if necessary and places implants
- Dentist places impression coping and takes an impression inside the immediate denture
- Lab converts immediate denture to Screw Retained Provisional
- Dentist screws Provisional in and adjusts occlusion



### ***WHAT DOES THE LAB NEED IN YOUR OFFICE?***

#### **WHAT DOES THE LAB NEED IN YOUR OFFICE?**

- 6 feet of countertop
- Sink with running water
- Electrical outlet
- Trash can



### ***WHAT DO YOU NEED FOR THE CONVERSION?***

#### **WHAT DO YOU NEED FOR THE CONVERSION?**

- Typical dental equipment
- Polyvinyl siloxane or polyether impression material
- Impression gun and tips
- Something to adjust occlusion
- Material to fill screw holes
- The Robertson Dental Lab Technician brings everything else