



209 W. Central Avenue • Lompoc, CA 93436
 805.735.7888 • 800.585.3111 • FAX: 805.735.2388
 Email: info@robertsondentallab.com
 www.robertsondentallab.com

WE ACCEPT DIGITAL IMPRESSIONS
 for design and/or complete fabrication.

Send STL files to
 INFO@RobertsonDentalLab.com

Call us or visit our website (resource page)
 for more information.

robertsondentallab.com

MATERIAL

Zirconia

- Monolithic MPa 1200
- Aesthetic-Ultra MPa 575
- Aesthetic-Super MPa 748
- Layered
 - Zirconia Occlusal/Lingual

Lithium Disilicate

- Monolithic Crown
- Layered Veneer

PFM (select metal below)

- FullCast (all metal units)**
 - Non-precious
 - Noble (semi-precious)*
 - High Noble White (precious)
 - High Noble Yellow (precious)

- Metal Try-in** Diagnostic Waxup Temporaries

DESIGN

Occlusion

- 0.00 mm
- 0.25 mm out*
- 0.50 mm out
- 1.00 mm out

Buccal Margin Design

- Supported Margin
- Porcelain Butt Margin


If No Occlusal Clearance

- Call Dr.
- Spot Prep
- Spot Opposing
- Metal Occlusal
- Make This My Preference
- Yes No

METAL DESIGN

- Coping With Full Porcelain Coverage
- Metal Coping With Porcelain Coverage*
- Coping 360° Collar 0.5mm
- Metal Occlusal Excluding Buccal Cusp
- Metal Occlusal Including Buccal Cusp
- Metal Lingual

Pontic Design



Occlusal Staining

- None
- Light*
- Medium
- Heavy



Stamp: **Shade:** **Finish:** **Try-in:** **Mould:**


Alameter Measurement _____

Papillameter Measurement _____

PLEASE INDICATE

Tooth Map Upper

Tooth Map Lower



*Standard unless specified otherwise.

Removable Prosthodontics

Preferred Tooth:

- Economy GENISYS Digital Denture Add AED
- Mid Denture
- Premium Immediate Denture

Partial Denture:

- Metal Frame Bite Rims
- Duraflex Surgical Guide
- Zirlux Acetal Custom Tray
- Stayplate/Flipper Mouthguard
- Wire Clasp Tap III
- Space Maintainer Repair/Reline/Rebase
- Valplast Other - Specify in Rx
- Snap-On Smile Bleaching Tray
- Antiseptic Tray
- To Completion

Night Guards:

- Thermo Form Processed Key Soft
- Hard Hard/Soft Hard Only

Terms: Net 30 with a service charge of 2% per month on charges over 30 days after statement date. Only if signed, construct and deliver the herein described dental restoration. Client shall be responsible for all fees, costs, charges, and expenses expended or incurred in any suit or action for collection of past due amounts or enforcement of provisions of this Agreement.

Dr. Name _____

Date Prepped _____

Date Needed Back _____ Time _____

Patient Name _____

Patient Info: M ___ F ___ Age _____

Tooth Number _____



Call me to discuss case

Dr. Signature _____

Lic. # _____

IMPLANTS

Brand / System _____

Size _____

Cement retained Screw retained/Hybrid

Implant Temp

Abutment Material Type

Titanium abutment Zirconia abutment

Anodize metal abutment

PLEASE SEND MORE

RX's Lab Boxes Shipping Labels

Photos

Memory Card / Thumb Drive / CD Enclosed

Prints Enclosed

E-mailed to photos@robertsondentallab.com

Please check All that Apply

Pre-Prep Model Included Opposing Model

Case Has Been Disinfected Doctor to Trim Die

Bite Registration Enclosed Study Models

Please include a study model on all work involving anterior teeth.


Questions, Please call 805-735-7888

Restoration Checklist and Helpful Suggestions for Optimum Restorations

- Send a bite registration with any cases involving the terminal tooth of the arch or where there is no posterior stop.
- Send a pre-op study model with all anterior cases when possible
- Check for adjacent teeth for proper contact areas, and adjust if necessary.
- Use a single, full arch tray for cases with 3 units or more.
- Take a shade before prepping teeth
- Securely wrap impressions and models, and enclose any restorations in a small box

PREPARATION GUIDELINES


VENEER PREPARATION



1.5 mm - 2.0 mm




STANDARD PORCELAIN LABIAL MARGIN PREPARATION




ONLAY

1.0 mm Wall Thickness



1.5 mm - 2.0 mm



ALL-CERAMIC/COMPOSITE PREPARATIONS

BUCCAL
LINGUAL
INTERPROXIMAL } 1.0 - 1.5 mm
INCISAL 1.5 - 2.0 mm



LABIAL
LINGUAL
INTERPROXIMAL } 1.0 - 1.5 mm
OCCLUSAL 1.5 - 2.0 mm




INLAY

2.0 mm min.



2.0 mm min.



No Cavo-Surface Bevels

Laboratory Schedule

Local Courier Service: We offer twice-a-day courier service between the cities of Santa Maria and Santa Barbara. Once-a-day we visit Ventura, Ojai, Oxnard, Santa Paula, Paso Robles, Atascadero, San Luis Obispo, Morro Bay, and Cayucos for case pickup and delivery. Call (805) 735-7888 for case pickup or more info.

Mail Service: We use UPS/FedEx to pickup and deliver cases statewide & nationally. For shipping supplies, call 800.585.3111.

Please note that the schedule listed below includes the number of days the case will be in the laboratory. Shipping days should be added onto lab days when determining when the case will get back to your office.

Denture Repairs & Reline:

Santa Maria through Santa Barbara only:

- Same Day Delivery: Case must be picked up by morning delivery driver to have it delivered back same day.
- Next Day Delivery: Case picked up by afternoon delivery driver will be delivered the next Morning by noon (Friday afternoon pickup, delivered Monday by noon).

Nipomo to Paso Robles / Ventura, Camarillo, Oxnard, Santa Paula only:

- Next Day Delivery: Case picked up by our delivery driver will be delivered early next Afternoon.

Crown & Bridge

Zirconia, 3 units or less.....	3 days
Lithium Disilicate, 3 units or less.....	3 days
Zirconia, 4 units or more.....	5 days
Lithium Disilicate, 4 units or more.....	5 days
PFM.....	10 days
Full Cast.....	10 days
Diagnostic Waxups.....	5 days
Temporary Crowns.....	3 days

Implants

Complete Implant.....	12 days
Implant Crown on Solid Abutment.....	10 days
Implant Crown on Solid Abutment.....	10 days
Implant Temp with Titanium Abutment.....	2 days

Removables

Partial Frame.....	10 days
Set Partial Denture.....	5 days
Process Partial Denture.....	5 days
Set and Wax Full Denture.....	5 days
Process Full Denture.....	5 days
DuraFlex, Visiclear, <i>includes try-in & finish</i>	10 days
Custom Tray.....	3 days
Occlusal Rims and Base Plates.....	3 days
Stayplate - 3 Teeth or Less.....	3 days
Stayplate - 4 Teeth or More.....	5 days
Valplast Process.....	5 days
Emergency Replacement Denture.....	5 days
Reset.....	5 days
Rebase.....	3 days

Guards/Splints/ Ortho

Guards.....	5 days
Splints.....	5 days
Sleep Apnea Appliance.....	10 days
Space Maintainer.....	10 days
Hawley Retainer.....	10 days

Aligners

Aligners, Digital Treatment Plan.....	3 days
Aligners, Aligner Fabrication.....	7 days

Snap-On Smile

Snap-On Smile.....	6 days
--------------------	--------

RUSH CASES!

WE CAN ACCOMMODATE RUSH CASES. CALL FOR TIMELINE AND RUSH FEES!

800.585.3111

Robertson Dental Lab is committed to providing you and your patients with high-quality, on-time dental restorations that meet your expectations and specifications every time. Our CDT's are here to help guide you through even the most difficult cases. Experience the Robertson Advantage today!